

SMALL PARTS KITS DIRECTIONS FOR USE

Refer to the driving console Operator's Manual for Small Parts Kit compatibility.

CAUTION: The Small Parts Kit Directions-for-Use are not intended to substitute for the necessity of reading and understanding the driving console Operator's Manual. The Operator's Manual, which is provided with the instrument, includes in-depth material intended to familiarize the Operating Room Staff with the controls and functions of the instrument.

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.

CAUTION: The equipment used in conjunction with the Alcon Small Parts Kit constitutes a complete surgical system. Use of disposables other than those of Alcon may affect system performance and create potential hazards, and if it is determined to have contributed to the malfunction of equipment under contract, could result in the avoidance of the contract and/or invoicing at prevailing hourly rates.

DESCRIPTION: Each Small Parts Kit contains sterile disposable phacoemulsification procedure supplies. The Small Parts Kit contains (1) infusion sleeve with BSI, (1) infusion sleeve, (1) I/A tip wrench, and (1) test chamber
OR

The Small Parts Kit contains (2) infusion sleeves without BSI, (1) I/A tip wrench, (1) test chamber, and may contain (1) phaco tip and wrench.

DIRECTION FOR USE: Refer to driving console Operator's Manual and procedural pack Directions for Use for setup instructions. Refer to Table below for recommended infusion sleeve and tip combinations, and incision sizes.

Infusion Sleeves/Type	Sleeve Color	Recommended Tips	Recommended Incision Size
1.1 mm High Infusion Sleeve	Light Blue	1.1 mm and 1.1 mm ABS [®] Tip, 1.1 mm Flared ABS [®] Tip, 1.1 mm Mackool [®] Flared ABS [®] Tip, 0.9 mm ALCON [®] UltraChopper [®] Tip 1.1 mm and 1.1 MI AquaLase [®] Tip Standard I/A, Silicone I/A, INTREPID [®] I/A Tip	3.2 mm
1.1 mm Sleeve	Dark Blue	1.1 mm and 1.1 mm ABS [®] Tip, 1.1 mm Flared ABS [®] Tip, 1.1 mm Mackool [®] Flared ABS [®] Tip, 0.9 mm ALCON [®] UltraChopper [®] Tip 1.1 mm and 1.1 MI AquaLase [®] Tip Standard I/A, Silicone I/A, INTREPID [®] I/A Tip	3.0 mm
1.1 mm Micro Sleeve	Blue/Green	1.1 mm and 1.1 mm ABS [®] Tip, 1.1 mm Flared ABS [®] Tip, 1.1 mm Mackool [®] Flared ABS [®] Tip, 0.9 mm ALCON [®] UltraChopper [®] Tip 1.1 MI AquaLase [®] Tip Standard I/A, Silicone I/A, INTREPID [®] I/A Tip	2.75 mm
1.1 mm Ultra Sleeve	Green	1.1 mm and 1.1 mm ABS [®] Tip, 1.1 mm Flared ABS [®] Tip, 1.1 mm Mackool [®] Flared ABS [®] Tip, 0.9 mm ALCON [®] UltraChopper [®] Tip 1.1 MI AquaLase [®] Tip Silicone I/A, INTREPID [®] I/A Tip	2.2 mm
1.1 mm Nano Sleeve	Light Green	1.1 mm Flared ABS [®] Tip, 0.9 mm ALCON [®] UltraChopper [®] Tip Silicone I/A, INTREPID [®] I/A Tip	1.8 mm
0.9 mm High Infusion Sleeve	Light Purple	0.9 mm and 0.9 mm ABS [®] MicroTip, 0.9 mm Tapered ABS [®] MicroTip, 0.9 mm Flared ABS [®] MicroTip, 0.9 mm Mackool [®] Flared ABS [®] MicroTip, 0.9mm ABS [®] Mini and Mini-Flared ABS [®] Tip 0.9 mm ALCON [®] UltraChopper [®] Tip Standard I/A, Silicone I/A, INTREPID [®] I/A Tip	3.2 mm
0.9 mm Micro Sleeve	Dark Purple	0.9 mm and 0.9 mm ABS [®] MicroTip, 0.9 mm Tapered ABS [®] MicroTip, 0.9 mm Flared ABS [®] MicroTip, 0.9 mm Mackool [®] Flared ABS [®] MicroTip, 0.9mm ABS [®] Mini and Mini-Flared ABS [®] Tip 0.9 mm ALCON [®] UltraChopper [®] Tip Standard I/A, Silicone I/A, INTREPID [®] I/A Tip	2.75 mm
0.9 mm Ultra Sleeve	Red	0.9 mm Flared ABS [®] MicroTip, 0.9mm ABS [®] Mini and Mini-Flared ABS [®] Tip 0.9 mm ALCON [®] UltraChopper [®] Tip 0.9mm ABS [®] INTREPID [®] Balanced Tip Silicone I/A, INTREPID [®] I/A Tip	2.2 mm
0.9 mm Nano Sleeve	Orange	0.9 mm Flared ABS [®] MicroTip, 0.9mm ABS [®] Mini and Mini-Flared ABS [®] Tip 0.9 mm ALCON [®] UltraChopper [®] Tip 0.9mm ABS [®] INTREPID [®] Balanced Tip Silicone I/A, INTREPID [®] I/A Tip	1.8 mm

WARNINGS!

- Improper usage or assembly could result in a potentially hazardous condition for the patient:
 - Mismatch of consumable components and/or use of settings not specially adjusted for a particular combination of consumable components may create potentially hazardous fluidics imbalance.
 - Use of non-approved handpieces may create potentially hazardous fluidics imbalance.
 - Use of incisions that are smaller than recommended can lead to mechanical and/or thermal damage to the eye tissues.
- If the handpiece test chamber is collapsed after tuning, there is a potential of low irrigation flow through the handpiece and may result in fluidic imbalance. This, in turn, may cause a shallowing or collapsing of the anterior chamber.
- Use of this product may require setting adjustments especially with the Ultra and Nano Sleeves. It is recommended that the Ultra and Nano Sleeves only be used with the INFINITI® Vision System, LAUREATE® World Phaco System, CONSTELLATION® Vision System, and CENTURION® Vision System as the fluidics performance of other consoles may be insufficient to perform Micro Co-axial cataract procedures. Prior to initial use, contact your Alcon Sales representative for in-service information. (Within the U.S. call 800-TO-ALCON or 817-293-0450. Outside of the U.S., contact your local Alcon Sales Representative.)
- Good clinical practice dictates the testing for adequate irrigation and aspiration flow prior to entering the eye.
- If any of the items in the Small Parts Kit is received in a defective condition, Alcon is to be notified immediately. Do not use any of the contents if the sterile package is damaged or the seal is broken in any way. In these cases, please contact:

By Phone: In USA (800) 757-9780
 Ask for Medical Safety
 International (817) 293-0450
 Or contact local Alcon Representative

By Mail: Alcon Research, Ltd
 Attention: Medical Safety (AB2-6)
 6201 South Freeway
 Fort Worth, TX 76134-2099
 USA

By E-mail: MedicalSafetyHouston@alconlabs.com

Each Small Parts Kit is identified by a lot number which provides traceability and should be given to Medical Safety Department when discussing the Small Parts Kit.

- Sterile disposable medical devices should not be re-used. Alcon assumes no responsibility for complications which may arise as a result of reuse or improper usage of any part of the Small Parts Kit.
 - Potential risk from reuse or reprocessing include: fluid path leaks or obstruction resulting in reduced fluidics performance, and foreign particle introduction into the eye.**

Definitions for symbols that may appear on product labels:

	SEE DIRECTIONS FOR USE		DOES NOT CONTAIN LATEX OR DRY NATURAL RUBBER		BATCH CODE
	USE BY : YEAR-MONTH	One of the following sterilization symbols will apply to this package:			MANUFACTURER
	SINGLE USE - DO NOT REUSE				STERILE - STERILIZED BY IRRADIATION
	DO NOT USE IF PACKAGE IS DAMAGED		STERILE - STERILIZED BY ETHYLENE OXIDE		CATALOG NUMBER
	AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY				PEEL POINT
	CAUTION: U.S. FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN				DO NOT RESTERILIZE

Alcon Laboratories (U.K.) Ltd.
 Frimley Business Park, Frimley, Camberley
 Surrey, GU16 7SR, United Kingdom

0123

Alcon[®]
 ALCON LABORATORIES, INC.
 6201 SOUTH FREEWAY
 FORT WORTH, TX 76134-2099 USA
 a Novartis company
 MADE IN USA

*Mackool is a trademark of
 Richard J. Mackool, M. D.

© 2008-2011, 2013 Novartis

316-2764-001 REV. R