


GRIESHABER®

DSP Aspheric Macular Lens

Now Available!



- PRECISE** - precision crafted
- COMPLIANT** - sterile for every procedure
- RELIABLE** - high quality
- CONVENIENT** - ready to use



- *Designed for improved visualization with enhanced contrast and crisp viewing throughout field of view*
- *Excellent self retaining qualities with more than 200mm² corneal surface contact area*
- *Six tabs stabilize but do not interfere with cannulas or globe rotation.*
- *Aspheric geometry offers improvements over traditional lens geometry to compensate for off-axis stereo viewing, defocus, loss of contrast and edge sharpness*



GRIESHABER® DSP Aspheric Macular Lens
618.60 6/box ~30° FOV -59D



Lens is delivered in sterile single use packaging

Contact your local representative at 1-800-TO-ALCON (1-800-862-5266)
alconretina.com

GRIESHABER REVOLUTION® DSP (6/box)



Style	20G	23G	25G	25+™
GRIESHABER® MAXGrip™ forceps	707.13	706.13	705.13	705.13P
End-grasping forceps	707.12	706.43	705.43	705.43P
ILM forceps	707.44	706.44	705.44	705.44P
Asymmetrical forceps		706.45	705.45	705.45P
Serrated forceps	707.47	706.47	705.47	705.47P
Curved scissors		706.52	705.52	705.52P
Vertical scissors	707.25	706.26		
Fiber optic forceps	708.12			

Reusable Handles

RENAISSANCE® Advanced
Designed for use with DSP tips only.
713.00.38 Red



NEW

GRIESHABER REVOLUTION®

712.00.41 Red
712.00.42 Blue



Sutherland NG

712.00.11 Red
712.00.12 Blue



RENAISSANCE® NG

712.00.31 Red
712.00.32 Blue




GRIESHABER® Advanced DSP Tips (6/box)



Style	20G	23G	25G	25+™
GRIESHABER® MAXGrip™ forceps	720.13	723.13	725.13	725.13P
End-grasping forceps	720.43	723.43	725.43	725.43P
ILM forceps	720.44	723.44	725.44	725.44P
Asymmetrical forceps		723.45	725.45	725.45P
Serrated forceps	720.47	723.47	725.47	725.47P
Curved scissors	720.52	723.52	725.52	725.52P
Vertical Scissors	720.26	723.26		

Advanced Backflush DSP—Passive or Active Aspirating (6/box)



Style	25G	23G	20G
Blunt	337.83	337.85	337.81
Soft Tip	337.84	337.86	337.82



Soft-Tip Cannulas (10/box)

Style	25G	23G	20G
Soft Tip with silicone end	8065149525	8065149523	8065149520

MVR Blade (6/box)



Style	25G	23G	20G
MVR Blade	8065912501	8065912301	8065912001

GRIESHABER® DSP Aspheric Macular Lens

618.60 6/box ~30° FOV -59D



Quick Lock Adaptor (reusable) 712.00.01 1 each

Converts the ACCURUS® Pneumatic Handpiece (8065808101) to be used with all DSP tips.



Products and packaging are color coded for easy identification.



Indications for Use: The GRIESHABER® DSP Aspheric Macular Lens is used to visualize the fundus and retinal structures during vitreoretinal surgery. It is designed as a self-retaining contact lens to allow hands-free operation.

Caution: Federal (USA) law restricts this device to sale by, or on the order of, a physician.

Warnings and Precautions:

- The device may be used only by well trained physicians and personnel.

- Potential risk from reuse or reprocessing include: reduced optical quality, surface damage on the optics, and foreign particle introduction to the eye.
- The surgeon's team has to ensure that sufficient viscous fluid is available throughout usage of the lens for continuous humidification of the cornea.

Attention: Reference the Directions for Use for a complete listing of indications, warnings, and precautions.