

AIR OPTIX® BRAND CONTACT LENSES US LENS PARAMETER GUIDE

November 2016



AIR OPTIX® plus HydraGlyde Contact Lenses

Base Curve(s):	Diameter(s):	Power Range(s):	Color(s):
8.6	14.2	+6.00 to -8.00 <i>(in 0.25 steps)</i> +6.50 to +8.00 -8.50 to -12.00 <i>(in 0.50 steps)</i>	Light Blue VISITINT®

Material: Silicone Hydrogel, lotrafilcon B **Water Content:** 33% **Dk/t:** 138 @-3.00D **Surface:** SmartShield® Technology
Center thickness: 0.08mm @-3.00D **Wearing Schedule:** Daily wear and up to 6 nights extended wear
Recommended replacement schedule: Monthly **Recommended lens care:** OPTI-FREE® PureMoist® Multi-Purpose Disinfecting Solution or CLEAR CARE® PLUS Cleaning & Disinfection Solution



AIR OPTIX® AQUA Contact Lenses

Base Curve(s):	Diameter(s):	Power Range(s):	Color(s):
8.6	14.2	+6.00 to -8.00 <i>(in 0.25 steps)</i> -8.50 to -10.00 <i>(in 0.50 steps)</i>	Light Blue VISITINT®

Material: Silicone Hydrogel, lotrafilcon B **Water Content:** 33% **Dk/t:** 138 @-3.00D **Surface:** SmartShield® Technology
Center thickness: 0.08mm @-3.00D **Wearing Schedule:** Daily wear and up to 6 nights extended wear
Recommended replacement schedule: Monthly **Recommended lens care:** OPTI-FREE® PureMoist® Multi-Purpose Disinfecting Solution or CLEAR CARE® Cleaning & Disinfection Solution



AIR OPTIX® NIGHT & DAY® AQUA Contact Lenses

Base Curve(s):	Diameter(s):	Power Range(s):	Color(s):
8.4, 8.6	13.8	+6.00 to -8.00 <i>(in 0.25 steps)</i> -8.50 to -10.00 <i>(in 0.50 steps)</i>	Light Blue VISITINT®

Material: Silicone Hydrogel, lotrafilcon A **Water Content:** 24% **Dk/t:** 175@-3.00D **Surface:** SmartShield® Technology
Center thickness: 0.08mm@-3.00D **Wearing schedule:** FDA-approved for daily wear and up to 30 nights of continuous wear
Recommended replacement schedule: Monthly **Indication:** Includes vision correction and therapeutic use
Recommended lens care: OPTI-FREE® PureMoist® Multi-Purpose Disinfecting Solution or CLEAR CARE® Cleaning & Disinfection Solution

PERFORMANCE DRIVEN BY SCIENCE®

AIR OPTIX® BRAND CONTACT LENSES US LENS PARAMETER GUIDE

November 2016



AIR OPTIX® COLORS Contact Lenses

Base Curve(s):	Diameter(s):	Power Range(s):	Color(s):
8.6	14.2	+6.00 to -6.00 (in 0.25 steps) -6.50 to -8.00 (in 0.50 steps)	Subtle Colors: Blue Green Brown Pure Hazel Gray Vibrant Colors: Brilliant Blue Gemstone Green Sterling Gray Honey

Material: Silicone Hydrogel, lotrafilcon B **Water Content:** 33% **Dk/t:** 138@-3.00D **Surface:** SmartShield® Technology
Center thickness: 0.08 mm@-3.00D **Wearing schedule:** Daily wear only **Recommended replacement schedule:** Monthly
Recommended lens care: OPTI-FREE® PureMoist® Multi-Purpose Disinfecting Solution or CLEAR CARE® Cleaning & Disinfection Solution



AIR OPTIX® for Astigmatism Contact Lenses

Base Curve(s):	Diameter(s):	Power Range(s):	Cylinder(s):	Axes:	Color(s):
8.7	14.5	+6.00 to -6.00 (in 0.25 steps) -6.50 to -10.00 (in 0.50 steps)	-0.75, -1.25, -1.75, -2.25	10 to 180 (in 10 degrees)	Light Blue VISITINT®

Material: Silicone Hydrogel, lotrafilcon B **Water Content:** 33% **Dk/t:** 108@-3.00 -125 x 180 **Surface:** SmartShield® Technology
Center thickness: 0.102mm@-3.00D **Design:** Back Surface Toric **Stabilization:** Unique precision balance 8 | 4™ design. Scribe marks at 3, 6, and 9 o'clock (6 o'clock scribe mark is slightly wider to differentiate from 3 and 9). **Wearing schedule:** Daily wear and up to 6 nights extended wear
Recommended replacement schedule: Monthly **Recommended lens care:** OPTI-FREE® PureMoist® Multi-Purpose Disinfecting Solution or CLEAR CARE® Cleaning & Disinfection Solution



AIR OPTIX® AQUA Multifocal Contact Lenses

Base Curve(s):	Diameter(s):	Power Range(s):	ADDs:	Color(s):
8.6	14.2	+6.00 to -10.00 (in 0.25 steps)	LO MED HI	Light Blue VISITINT®

Material: Silicone Hydrogel, lotrafilcon B **Water Content:** 33% **Dk/t:** 138@-3.00D **Surface:** SmartShield® Technology
Center thickness: 0.08mm@-3.00D **Wearing schedule:** Daily wear and up to 6 nights extended wear
Recommended replacement schedule: Monthly **Recommended lens care:** OPTI-FREE® PureMoist® Multi-Purpose Disinfecting Solution or CLEAR CARE® Cleaning & Disinfection Solution

PERFORMANCE DRIVEN BY SCIENCE®

DAILIES® BRAND CONTACT LENSES US LENS PARAMETER GUIDE

November 2016



DAILIES® TOTAL1® Contact Lenses

Base Curve(s):	Diameter(s):	Power Range(s):	Color(s):
8.5	14.1	+0.50 to +6.00 -0.50 to -6.00 <i>(in 0.25 steps)</i> -6.50 to to -12.00 <i>(in 0.50 steps)</i>	Light Blue VISITINT®

Material: Water Gradient, delefilcon A **Water Content:** 33% **Surface Water Content:** ≥80% **Dk/t:** 156@ -3.00
Center thickness: 0.09 mm@-3.00D (varies with power) **Wearing schedule:** Daily wear **Recommended replacement schedule:** Daily disposable
Recommended lens care: The lenses are not intended to be cleaned or disinfected and should be discarded after a single use.



DAILIES® TOTAL1® Multifocal Contact Lenses

Base Curve(s):	Diameter(s):	Power Range(s):	Color(s):
8.5	14.1	+6.00 to -10.00 <i>(in 0.25 steps)</i>	Light Blue VISITINT®

Material: Water Gradient, delefilcon A **Water Content:** 33% **Surface Water Content:** ≥80% **Dk/t:** 156@ -3.00
Center thickness: 0.09 mm@-3.00D (varies with power) **Wearing schedule:** Daily wear **Recommended replacement schedule:** Daily disposable
Recommended lens care: The lenses are not intended to be cleaned or disinfected and should be discarded after a single use.



DAILIES® AquaComfort Plus® Contact Lenses

Base Curve(s):	Diameter(s):	Power Range(s):	Color(s):
8.7	14.0	+0.50 to +6.00 -0.50 to -6.00 <i>(in 0.25 steps)</i> +6.50 to +8.00 -6.50 to -15.00 <i>(in 0.50 steps)</i>	Light Blue VISITINT®

Material: nelfilcon A. **Water Content:** 69% **Dk/t:** 26@ -3.00 **Center thickness:** 0.10mm@-3.00D (varies with power) **Wearing schedule:** Daily wear
Recommended replacement schedule: Daily disposable **Recommended lens care:** The lenses are not intended to be cleaned or disinfected and should be discarded after a single use.



DAILIES® AquaComfort Plus® Toric Contact Lenses

Base Curve(s):	Diameter(s):	Power Range(s):	Cylinder(s):	Axes:	Color(s):
8.8	14.4	Plano to -6.00 +0.25 to +4.00 <i>(in 0.25 steps)</i> -6.50 to -8.00 <i>(in 0.50 steps)</i>	-0.75, -1.25, -1.75	10/20/70/80/ 90/100/110/ 160/170/180	Light Blue VISITINT®

Material: nelfilcon A. **Water Content:** 69% **Dk/t:** 26@ -3.00 **Center thickness:** 0.10mm@-3.00D (varies with power)
Wearing schedule: Daily wear **Recommended replacement schedule:** Daily disposable **Recommended lens care:** The lenses are not intended to be cleaned or disinfected and should be discarded after a single use.



DAILIES® AquaComfort Plus® Multifocal Contact Lenses

Base Curve(s):	Diameter(s):	Power Range(s):	ADDs:	Color(s):
8.7	14.0	Plano to -10.00 +0.50 to +6.00 <i>(in 0.25 steps)</i>	LO MED HI	Light Blue VISITINT®

Material: nelfilcon A. **Water Content:** 69% **Dk/t:** 26@ -3.00 **Center thickness:** 0.10mm@-3.00D (varies with power) **Wearing schedule:** Daily wear
Recommended replacement schedule: Daily disposable **Recommended lens care:** The lenses are not intended to be cleaned or disinfected and should be discarded after a single use.

PERFORMANCE DRIVEN BY SCIENCE®

FreshLook® BRAND CONTACT LENSES US LENS PARAMETER GUIDE

November 2016



FRESHLOOK® HANDLING TINT UV Contact Lenses

Base Curve(s):	Diameter(s):	Power Range(s):	Color(s):
Median	14.5	-0.25 to -6.00 +0.25 to +2.00 <i>(in 0.25 steps)</i> -6.50 to -8.00 +2.50 to +6.00 <i>(in 0.50 steps)</i>	Light Green Handling Tint

Material: phemfilcon A. **Water Content:** 55% **Dk/l:** 20@ -3.00 **Center thickness:** 0.08mm **Wearing schedule:** Daily wear only
Recommended replacement: 2 weekly **Recommended lens care:** OPTI-FREE® PureMoist® Multi-Purpose Disinfecting Solution or CLEAR CARE® Cleaning & Disinfection Solution



FreshLook® COLORBLEND® Contact Lenses

Base Curve(s):	Diameter(s):	Power Range(s):	Color(s):
Median	14.5	Plano to -6.00 +0.25 to +2.00 <i>(in 0.25 steps)</i> -6.50 to -8.00 +2.50 to +6.00 <i>(in 0.50 steps)</i>	Amethyst Blue Brilliant Blue Brown Gemstone Green Gray Green Honey Pure Hazel Sterling Gray True Sapphire Turquoise

Material: phemfilcon A. **Water Content:** 55% **Dk/l:** 20@ -3.00 **Center thickness:** 0.08mm **Wearing schedule:** Daily wear only
Recommended replacement schedule: 2 weekly **Recommended lens care:** OPTI-FREE® PureMoist® Multi-Purpose Disinfecting Solution or CLEAR CARE® Cleaning & Disinfection Solution



FreshLook® COLORS Contact Lenses

Base Curve(s):	Diameter(s):	Power Range(s):	Color(s):
Median	14.5	Plano to -6.00 +0.25 to +2.00 <i>(in 0.25 steps)</i> -6.50 to -8.00 +2.50 to +6.00 <i>(in 0.50 steps)</i>	Blue Green Hazel Violet

Material: phemfilcon A. **Water Content:** 55% **Dk/l:** 20@ -3.00 **Center thickness:** 0.08mm **Wearing schedule:** Daily wear only
Recommended replacement schedule: 2 weekly **Recommended lens care:** OPTI-FREE® PureMoist® Multi-Purpose Disinfecting Solution or CLEAR CARE® Cleaning & Disinfection Solution



FreshLook® DIMENSIONS Contact Lenses

Base Curve(s):	Diameter(s):	Power Range(s):	Color(s):
Median	14.5	Plano to -6.00 +0.25 to +2.00 <i>(in 0.25 steps)</i> -6.50 to -8.00 +2.50 to +6.00 <i>(in 0.50 steps)</i>	Caribbean Aqua Pacific Blue Sea Green

Material: phemfilcon A. **Water Content:** 55% **Dk/l:** 20@ -3.00 **Center thickness:** 0.08mm **Wearing schedule:** Daily wear only
Recommended replacement schedule: 2 weekly **Recommended lens care:** OPTI-FREE® PureMoist® Multi-Purpose Disinfecting Solution or CLEAR CARE® Cleaning & Disinfection Solution



FreshLook® One-Day Contact Lenses

Base Curve(s):	Diameter(s):	Power Range(s):	Color(s):
8.6	13.8	Plano, -0.50 to -6.00 <i>(in 0.25 steps)</i> -6.50 to -8.00 <i>(in 0.50 steps)</i>	Blue Green Pure Hazel Gray

Material: nelfilcon A. **Water Content:** approximately 69% **Dk/l:** 26@ -3.00 **Center thickness:** 0.10mm@ -3.00
Wearing schedule: Daily disposable **Recommended replacement schedule:** Daily **Recommended lens care:** The lenses are not intended to be cleaned or disinfected and should be discarded after a single use.

PERFORMANCE DRIVEN BY SCIENCE®

FOCUS® DAILIES® BRAND CONTACT LENSES US LENS PARAMETER GUIDE

November 2016



FOCUS® DAILIES® Contact Lenses

Base Curve(s):	Diameter(s):	Power Range(s):	Color(s):
8.6	13.8	-0.50 to -6.00 <i>(in 0.25 steps)</i> -6.50 to -10.00 <i>(in 0.50 steps)</i> +0.50 to +6.00 <i>(in 0.25 steps)</i>	Light Blue VISITINT®

Material: nelfilcon A. **Water Content:** 69% **Dk/t:** 26 @ -3.00 **Center thickness:** 0.10 mm at -3.00D (varies with power)

Wearing schedule: Daily wear **Recommended replacement:** Daily Disposable.

Recommended lens care: The lenses are not intended to be cleaned or disinfected and should be discarded after a single use.



FOCUS® DAILIES® Toric Contact Lenses

Base Curve(s):	Diameter(s):	Power Range(s):	Cylinder(s):	Axes:	Color(s):
8.6	14.2	+4.00 to -6.00 <i>(in 0.25 steps)</i> -6.50 to -8.00 <i>(in 0.50 steps)</i>	-0.75, -1.50	20/70/90 110/160/180	Light Blue VISITINT®

Material: nelfilcon A. **Water Content:** 69% **Dk/t:** 26 @ -3.00 **Center thickness:** 0.10 mm at -3.00D (varies with power)

Wearing schedule: Daily wear **Recommended replacement:** Daily Disposable

Recommended lens care: The lenses are not intended to be cleaned or disinfected and should be discarded after a single use.

PERFORMANCE DRIVEN BY SCIENCE®

Important information for AIR OPTIX® plus HydraGlyde (lotrafilcon B) contact lenses: For daily wear or extended wear up to 6 nights for near/far-sightedness. Risk of serious eye problems (i.e., corneal ulcer) is greater for extended wear. In rare cases, loss of vision may result. Side effects like discomfort, mild burning or stinging may occur.

Important information for AIR OPTIX® AQUA (lotrafilcon B) contact lenses: For daily wear or extended wear up to 6 nights for near / far-sightedness. Risk of serious eye problems (i.e., corneal ulcer) is greater for extended wear. In rare cases, loss of vision may result. Side effects like discomfort, mild burning or stinging may occur.

Important information for AIR OPTIX® COLORS (lotrafilcon B) contact lenses: For daily wear only for near / far-sightedness. Contact lenses, even if worn for cosmetic reasons, are prescription medical devices that must only be worn under the prescription, direction and supervision of an eye care professional. Serious eye health problems may occur as a result of sharing contact lenses. Although rare, serious eye problems can develop while wearing contact lenses. Side effects like discomfort, mild burning or stinging may occur. To help avoid these problems, patients must follow the wear and replacement schedule and the lens care instructions provided by their eye doctor.

Important information for AIR OPTIX® for Astigmatism (lotrafilcon B) contact lenses: For daily wear or extended wear up to 6 nights for near / far-sightedness and astigmatism. Risk of serious eye problems (i.e., corneal ulcer) is greater for extended wear. In rare cases, loss of vision may result. Side effects like discomfort, mild burning or stinging may occur.

Important information for AIR OPTIX® AQUA Multifocal (lotrafilcon B) contact lenses: For daily wear or extended wear up to 6 nights for near / far-sightedness and/or presbyopia. Risk of serious eye problems (i.e., corneal ulcer) is greater for extended wear. In rare cases, loss of vision may result. Side effects like discomfort, mild burning or stinging may occur.

Important information for AIR OPTIX® NIGHT & DAY® AQUA (lotrafilcon A) contact lenses: Indicated for vision correction for daily wear (worn only while awake) or extended wear (worn while awake and asleep) for up to 30 nights. Relevant Warnings: A corneal ulcer may develop rapidly and cause eye pain, redness or blurry vision as it progresses. If left untreated, a scar, and in rare cases loss of vision, may result. The risk of serious problems is greater for extended wear vs. daily wear and smoking increases this risk. A one-year post-market study found 0.18% (18 out of 10,000) of wearers developed a severe corneal infection, with 0.04% (4 out of 10,000) of wearers experiencing a permanent reduction in vision by two or more rows of letters on an eye chart. Relevant Precautions: Not everyone can wear for 30 nights. Approximately 80% of wearers can wear the lenses for extended wear. About two-thirds of wearers achieve the full 30 nights continuous wear. Side Effects: In clinical trials, approximately 3-5% of wearers experience at least one episode of infiltrative keratitis, a localized inflammation of the cornea which may be accompanied by mild to severe pain and may require the use of antibiotic eye drops for up to one week. Other less serious side effects were conjunctivitis, lid irritation or lens discomfort including dryness, mild burning or stinging. Contraindications: Contact lenses should not be worn if you have: eye infection or inflammation (redness and/or swelling); eye disease, injury or dryness that interferes with contact lens wear; systemic disease that may be affected by or impact lens wear; certain allergic conditions or using certain medications (ex. some eye medications). Additional Information: Lenses should be replaced every month. If removed before then, lenses should be cleaned and disinfected before wearing again. Always follow the eye care professional's recommended lens wear, care and replacement schedule. Consult package insert for complete information, available without charge by calling (800) 241-5999 or go to myalcon.com

See product instructions for complete wear, care and safety information. 

© 2016 Novartis 10/16 US-DAF-15-E-0876(2)

Alcon A Novartis
Division